

Anti-CD137/4-1BB/TNFRSF9 antibody [3H3] {PE} (STJA0031458)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | PE |
| Purification | Purified |
| Formulation | PBS pH 7.4, 50% glycerol, 0.05% proclin300 |
| Isotype | IgG2ak |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 21942 |
| Gene Symbol | Tnfrsf9 |
| UniProt ID | TNR9_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|