

Anti-AMHR2/MISRll antibody [12G4] {PE} (STJA0031522)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/IF/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PE
Purification	Purified
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	269
Gene Symbol	AMHR2
UniProt ID	AMHR2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--