

Anti-CD344/FZD4 antibody [5017#] {PE} (STJA0031646)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Purification	Purified
Formulation	0.2M Borate buffer pH8.5, 0.02% BSA, 0.05% Proclin 300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	8322
Gene Symbol	FZD4
UniProt ID	FZD4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--