

**Anti-CD20/MS4A1 antibody [SAA2209] {FITC} (STJA0031790)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Human
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	FITC
Purification	Affinity chromatography.
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">931</a>
Gene Symbol	<a href="#">MS4A1</a>
UniProt ID	<a href="#">CD20_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------