

**Anti-CD19 antibody [SAA2206] {APC} (STJA0031811)**

| GENERAL INFORMATION    |                     |                                                                                                                      |
|------------------------|---------------------|----------------------------------------------------------------------------------------------------------------------|
|                        | Product Type        | Primary antibodies                                                                                                   |
|                        | Applications        | ELISA/FC                                                                                                             |
|                        | Host / Source       | Human                                                                                                                |
|                        | Reactivity          | Human                                                                                                                |
|                        |                     |                                                                                                                      |
| PRODUCT PROPERTIES     |                     |                                                                                                                      |
|                        | Clonality           | Monoclonal                                                                                                           |
|                        | Conjugation         | APC                                                                                                                  |
|                        | Purification        | Affinity chromatography.                                                                                             |
|                        | Formulation         | PBS pH 7.4, 50% glycerol, 0.05% proclin300                                                                           |
|                        | Isotype             | IgG1k                                                                                                                |
|                        | Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze.                                                       |
|                        |                     |                                                                                                                      |
| TARGET INFORMATION     |                     |                                                                                                                      |
|                        | Gene ID             | <a href="#">930</a>                                                                                                  |
|                        | Gene Symbol         | <a href="#">CD19</a>                                                                                                 |
|                        | UniProt ID          | <a href="#">CD19_HUMAN</a>                                                                                           |
|                        |                     |                                                                                                                      |
| ADDITIONAL INFORMATION |                     |                                                                                                                      |
|                        | Note                | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |