

**Anti-CD13/ANPEP/APN antibody [3L94] {PerCP} (STJA0031822)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/IF/WB
Host / Source	Mouse
Reactivity	Pig

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	PerCP
Purification	Affinity chromatography.
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	VHH-mFc (IgG2a)
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--