

Anti-CD227/MUC1 antibody [SAA2216] {FITC} (STJA0031834)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/FC
	Host / Source	Human
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	FITC
	Purification	Affinity chromatography.
	Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
	Isotype	IgG1k
	Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.
TARGET INFORMATION		
	Gene ID	4582
	Gene Symbol	MUC1
	UniProt ID	MUC1_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.