

## Anti-CD77 antibody [5B5] {PerCp} (STJA0031851)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PerCp
Purification	Affinity chromatography.
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--