

Anti-CD227/MUC1 antibody [HMFG1] {PerCp} (STJA0031857)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC/ICC/IF/IHC/WB
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PerCp
Purification	Protein A/G purified from cell culture supernatant.
Formulation	PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--