

**Anti-CD8A antibody [FS0307] {FITC} (STJA0031903)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FCM
Host / Source	Mouse
Reactivity	Bovine

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	FITC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--