

Anti-CD45/PTPRC antibody [BC8] {AF488} (STJA0031931)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FCM
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	AF488
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG1k
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
UniProt ID	PTPRC_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--