

Anti-CD185/CXCR5 antibody [Fs0373] {PE} (STJA0031974)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FCM |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------------------------------------|
| Clonality | Monoclonal |
| Conjugation | PE |
| Purification | Affinity chromatography |
| Formulation | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300. |
| Isotype | IgG1k |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 643 |
| Gene Symbol | CXCR5 |
| UniProt ID | CXCR5_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|