

## Anti-CD94/KLRD1 antibody [Fs0381] {PerCP} (STJA0031981)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FCM/FuncS
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PerCP
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">3824</a>
Gene Symbol	<a href="#">KLRD1</a>
UniProt ID	<a href="#">KLRD1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--