

**Anti-CD94/KLRD1 antibody [Fs0381] {AF488} (STJA0031982)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FCM/FuncS
Host / Source	Human
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	AF488
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG1k
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">3824</a>
Gene Symbol	<a href="#">KLRD1</a>
UniProt ID	<a href="#">KLRD1_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--