

## Anti-HLA-B\*27 antibody [HLA-ABC-m3] {FITC} (STJA0032004)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FCM
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	FITC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">3106</a>
Gene Symbol	<a href="#">HLA-B</a>
UniProt ID	<a href="#">HLAB_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--