

**Anti-CSF2/GM-CSF antibody [Fs0380] {AF488} (STJA0032014)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FCM/FuncS    |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

**PRODUCT PROPERTIES**

|                     |                                                                |
|---------------------|----------------------------------------------------------------|
| Clonality           | Monoclonal                                                     |
| Conjugation         | AF488                                                          |
| Purification        | Affinity chromatography                                        |
| Formulation         | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.                |
| Isotype             | IgG2ak                                                         |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">12981</a>      |
| Gene Symbol | <a href="#">Csf2</a>       |
| UniProt ID  | <a href="#">CSF2_MOUSE</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|