

Anti-EPOR antibody [Fs0378] {PerCP} (STJA0032019)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FCM
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PerCP
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	13857
Gene Symbol	Epor
UniProt ID	EPOR_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--