

Anti-IL-12 p40 antibody [Fs0379] {APC} (STJA0032023)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FCM |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | APC |
| Purification | Affinity chromatography |
| Formulation | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300. |
| Isotype | IgG2ak |
| Storage Instruction | Store at 4°C for 12 months. Protect from light. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 16160 |
| Gene Symbol | Il12b |
| UniProt ID | IL12B_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|