

Anti-IL-12 p40 antibody [Fs0379] {AF488} (STJA0032026)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FCM
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	AF488
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG2ak
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	16160
Gene Symbol	Il12b
UniProt ID	IL12B_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--