

Anti-CXCL9/MIG antibody [MIG-2F5.5] {FITC} (STJA0032028)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FCM/FuncS
Host / Source	Armenian hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	FITC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

TARGET INFORMATION

Gene ID	17329
Gene Symbol	Cxcl9
UniProt ID	CXCL9_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--