

**Anti-CXCL9/MIG antibody [MIG-2F5.5] {PerCP} (STJA0032031)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FCM/FuncS
Host / Source	Armenian hamster
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	PerCP
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">17329</a>
Gene Symbol	<a href="#">Cxcl9</a>
UniProt ID	<a href="#">CXCL9_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--