

Anti-CD3E antibody [Fs0316] {FITC} (STJA0032040)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FCM
Host / Source	Mouse
Reactivity	Cynomolgus monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	FITC
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months. Protect from light. Do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--