

## Anti-IA-2/PTPRN antibody [1A469] (STJA0032081)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IF/IHC/WB    |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Purification        | Protein A/G purified from cell culture supernatant   |
| Formulation         | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.  |
| Isotype             | IgG  |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">5798</a>        |
| Gene Symbol | <a href="#">PTPRN</a>       |
| UniProt ID  | <a href="#">PTPRN_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|