

Anti-GGGGS Linker antibody [R4E06] (STJA0032100)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FCM/ICC/IF/IHC/WB
Host / Source	Rabbit
Reactivity	General

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	PBS with 0.09% Sodium Azide
Isotype	IgG
Storage Instruction	Store at 2-8 °C 12 months.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--