

Anti-CKMB antibody [R4D88] (STJA0032106)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA-cap/CLIA/LF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

TARGET INFORMATION

Gene ID	1158 1152
Gene Symbol	CKM CKB
UniProt ID	KCRM_HUMAN KCRB_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--