

## Anti-Desmoglein 3 antibody [AK23] (STJA0032120)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IP
Host / Source	Mouse
Reactivity	Human/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2ak
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

### TARGET INFORMATION

Gene ID	<a href="#">1830</a>
Gene Symbol	<a href="#">DSG3</a>
UniProt ID	<a href="#">DSG3_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--