

## Anti-IL24 antibody [Fs0359] (STJA0032125)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FCM/Neut
Host / Source	Human
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG4k
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

### TARGET INFORMATION

Gene ID	<a href="#">11009</a>
Gene Symbol	<a href="#">IL24</a>
UniProt ID	<a href="#">IL24_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--