

Anti-ZEBOV GP/Envelope glycoprotein antibody [KZ52] (STJA0032150)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IF/Neut
Host / Source	Human
Reactivity	Zaire ebolavirus strain Mayinga-76/ZEBOV/Zaire Ebola virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--