

Anti-Coxsackievirus A16 Mature Virion/VP1, 2, 3 antibody [9D7] (STJA0032165)**GENERAL INFORMATION**

| | |
|---------------|--|
| Product Type | Primary antibodies |
| Applications | ELISA/Neut/XRC |
| Host / Source | Human |
| Reactivity | Coxsackievirus A16 strain Tainan 5079 98 |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | IgG1k |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|