

Anti-Flavivirus/YFV Envelope Protein E antibody [Z6] (STJA0032177)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FCM/Neut
Host / Source	Human
Reactivity	Yellow fever viruses/YFV/Zika virus/ZIKV/Dengue virus type 1/DENV-1/Dengue virus type 2/DENV-2/Dengue virus type 3/DENV-3/Dengue virus type 4/DENV-4

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -85°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--