

Anti-YFV Envelope protein E antibody [Fs0370] (STJA0032186)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/Neut
Host / Source	Mouse
Reactivity	DENV-1/Dengue virus type 1/DENV-2/Dengue virus type 2/DENV-3/Dengue virus type 3/DENV-4/Dengue virus type 4/ZIKV/Zika virus/JEV/Japanese encephalitis virus/Yellow Fever Virus/YFV

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2ak
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--