

Anti-HPIV3 HN/Hemagglutinin-neuraminidase antibody [Fs0339] (STJA0032188)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut
Host / Source	Human
Reactivity	Human parainfluenza 3 virus strain Wash 47885 57/HPIV-3/Human parainfluenza 3 virus strain NIH 47885

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--