

Anti-AAV9 VP1/Capsidprotein VP1 antibody (STJA0032195)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IFA/IHC/CLIA/LFA/IP
Host / Source	Mouse
Reactivity	Adeno-Associated Virus 9/AAV9

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--