

## Anti-Francisella tularensis LPS antibody [FB11] (STJA0032213)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/ICC/IF/WB
Host / Source	Mouse
Reactivity	F.tularensis

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300
Isotype	IgG2a
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

### TARGET INFORMATION

Specificity	There is no crossreactivity with Fr. novicida, Br. abortus, Br. suis, Br. melitensis, Br. ovis, Y. pestis, Y. enterocolitica, Y. pseudotuberculosis, E. coli, and V. cholerae.
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--