

Anti-CD20 & CD3e Bispecific antibody [Iv0225] (STJA0032269)

GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A or G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.

TARGET INFORMATION

Gene ID	12482 12501
Gene Symbol	Ms4a1 Cd3e
UniProt ID	CD20_MOUSE CD3E_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--