

## Anti-CD19 & CD3e Bispecific antibody [Iv0224] (STJA0032275)

### GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	Mouse
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A or G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.

### TARGET INFORMATION

Gene ID	<a href="#">12478</a> <a href="#">12501</a>
Gene Symbol	<a href="#">Cd19</a> <a href="#">Cd3e</a>
UniProt ID	<a href="#">CD19_MOUSE</a> <a href="#">CD3E_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--