

Anti-Lujo virus/LUJV Glycoprotein/GPC antibody [Fs0336] (STJA0032333)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut
Host / Source	Human
Reactivity	Lujo virus/LUJV

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol.
Isotype	VHH-hFc
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--