

Anti-SOX2 antibody [2A10-F11-F12] (STJA0032351)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IHC-IF/FC
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 ICC 1:150 FC 1:100
Formulation	Liquid in PBS with 0.09% (W/V) sodium azide
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	6657
Gene Symbol	SOX2
UniProt ID	SOX2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---