

Anti-COX IV antibody [4D11-B3-E8] (STJA0032352)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP/FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Hamster/Goat/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000 IHC 1:50-1:200 IF 1:100-1:200 IP 1:20 FC 1:100
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1327
Gene Symbol	COX4I1
UniProt ID	COX41_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




COX IV antibody western blot mouse skeletal muscle Jurkat STJA0032352




COX IV antibody immunocytochemistry HeLa cells STJA0032352




COX IV antibody immunoprecipitation HeLa lysates STJA0032352



COX IV antibody immunohistochemistry human colorectal carcinoma 4D11B3E8 STJA0032352



COX4 antibody flow cytometry K562 cells STJA0032352



COX IV antibody western blot COS7 3T3 K562 C6 CHO-K1 HeLa STJA0032352