

## Anti-PRMT6 antibody [2C3-F1-G1] (STJA0032353)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-F/IHC-P/ICC/IF |
| Host / Source | Mouse                 |
| Reactivity    | Human                 |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 1 mg/mL                                                                                  |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IF 1:50-1:200                                       |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">55170</a>      |
| Gene Symbol | <a href="#">PRMT6</a>      |
| UniProt ID  | <a href="#">ANM6_HUMAN</a> |

### ADDITIONAL INFORMATION


**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



PRMT6 antibody western blot HeLa 293T cells  
STJA0032353



PRMT6 antibody immunofluorescence HeLa cells  
STJA0032353



PRMT6 antibody immunohistochemistry human  
breast cancer STJA0032353