

Anti-SAFB1 antibody [2E8-E2-G6] (STJA0032393)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF/ChIP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200 ChIP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	6294
Gene Symbol	SAFB
UniProt ID	SAFB1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SAFB1 antibody western blot A431 MDA-MB-468 ZR751 MCF-7 STJA0032393



SAFB1 antibody immunocytochemistry HeLa cells STJA0032393