

Anti-PGP9.5 antibody [7B8-G5-B8] (STJA0032421)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Mouse
Reactivity	Human/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	3.1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	7345
Gene Symbol	UCHL1
UniProt ID	UCHL1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



UCHL1 PGP95 antibody western blot U251 UCHL1 PGP95 antibody immunocytochemistry
 IMR32 Raji COS7 STJA0032421 COS7 STJA0032421

