

## Anti-Hsp27 antibody [7E5-1F9-D9] (STJA0032423)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Mouse
Reactivity	Human/Monkey

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.2 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">3315</a>
Gene Symbol	<a href="#">HSPB1</a>
UniProt ID	<a href="#">HSPB1_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Hsp27 antibody western blot HeLa A431 COS7 LNCaP MCF-7 STJA0032423



Hsp27 antibody immunocytochemistry HeLa STJA0032423