

Anti-AMPK beta 1 antibody [1A7-E11-E9] (STJA0032433)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.47 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2a
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5564
Gene Symbol	PRKAB1
UniProt ID	AAKB1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




AMPK beta 1 antibody western blot multiple cell lines STJA0032433



AMPK beta 1 antibody immunocytochemistry HeLa STJA0032433



AMPK beta 1 antibody immunoprecipitation HeLa STJA0032433



AMPK beta 1 antibody immunohistochemistry breast cancer STJA0032433