

Anti-PI3 Kinase p85 beta antibody [8D9-D5-F8] (STJA0032438)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB
	Host / Source	Mouse
	Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5296
Gene Symbol	PIK3R2
UniProt ID	P85B_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



PI3 Kinase p85 beta antibody western blot
multiple cell lines STJA0032438