

## Anti-GFP/eGFP-Tag antibody [3F11-B12-E7-F12] (STJA0032442)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Species-independent

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.28 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



EGFP/EYFP antibody western blot Rosetta E coli STJA0032442



EGFP/EYFP antibody western blot CHO-K1 cells STJA0032442