

## Anti-STAT3 antibody [4C8-1C9-H8] (STJA0032449)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Monkey/Hamster

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	12.61 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6774</a>
Gene Symbol	<a href="#">STAT3</a>
UniProt ID	<a href="#">STAT3_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STAT3 antibody western blot CHOK1 MCF7 A549  
 C6 3T3 K562 HeLa STJA0032449