

Anti-Hsp90 beta antibody [3B9-D4-G4] (STJA0032457)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Monkey/Hamster

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	3326
Gene Symbol	HSP90AB1
UniProt ID	HS90B_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Hsp90 beta antibody western blot HeLa 3T3 C6 CHO-K1 COS7 STJA0032457



Hsp90 beta antibody immunocytochemistry HeLa STJA0032457



Hsp90 beta antibody immunoprecipitation HeLa lysates STJA0032457