

Anti-RSK1 p90 antibody [6B9-D6-F8] (STJA0032480)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	3.02 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	6195
Gene Symbol	RPS6KA1
UniProt ID	KS6A1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



p90RSK antibody western blot COS7 3T3 Jurkat HeLa A549 STJA0032480



p90RSK antibody immunoprecipitation HeLa lysates STJA0032480