

Anti-Lamin B1 antibody [5G8-D3-H7] (STJA0032486)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.46 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000 IF 1:100-1:200
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	4001
Gene Symbol	LMNB1
UniProt ID	LMNB1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Lamin B1 antibody western blot COS7 Molt4 A549 3T3 STJA0032486



Lamin B1 antibody immunocytochemistry HeLa cells STJA0032486